

[Wooju Mulsan Co.,Ltd]

Item	Steam Changing Eye Mask	Model	.
------	-------------------------	-------	---



Composition	Dimension(mm)	180*75 / 18g		
	Pcs/ Case	5pcs / Case		
	1Carton(pcs)	200	Carton Size(mm) / Kg	555*415*220 / 7kg

Description	<p><Characteristics></p> <ul style="list-style-type: none"> ■ Thermal interaction of Eye Mask helps your eye more comfortable. ■ You can recover skin of your eyes as the elastic twenties by using eye cream and eye mask at tired delicate skin around your eyes for a day. ■ The fumigation effect moisturizes and comforts your dry eyes. ■ Eye mask which is added ten kinds of Korean medicinal herbs makes the function and perfume better and human's body healthy. ■ Eye mask has soft wearable feeling to your eyes and makes eyes comfortable by using pure cotton cover . ■ Twenty Eye Cream contains Adenosine and Niasinamid ingredient to improve wrinkle and to help skin whitening and makes your eyes wet and elastic. <p><How to use></p> <ol style="list-style-type: none"> 1. Open eye cream pack, After cleaning face and arrange the skin with skin care . 2. Apply eye cream spread to your both around eyes from outside to inside as if massaging. 3. Open the Eye Mask and wearing (Thermal interaction is occurred on opening). 4. Take off Eye Mask 20 minutes later and check how much absorption of eye cream is occurred in your skin of around your eyes. The warm time of the eye mask lasts 2 to 3 hours. <p><Effect Function></p> <p>Eliminate dark eye circle, Relieve eye fatigue and strain, Help you sleep better & relieve dry eyes, Disposable mask makes it convenient and hygienic to use, Constant temperature heating, safe guarantee.</p> <p><Raw materials></p> <p>Heating element : Iron powder, activated carbon, water, inorganic salts etc. Pure cotton</p> <p><Notice></p> <ol style="list-style-type: none"> 1.Those who follow to below are prohibited before using. (The aged, infants, the sensitive to thermal. / under treatment from Doctor on account of eye disease / Those who have scars around the eyes or eczema etc. and other problem.) 2. Be careful to use high temp. heat elements, which cause cool burn without consequence. (Erythema , Blister) 3. Try not to use when sleeping for long time because of cold burn warning 4. Throw it in trash can after using. It is a disposable product and cannot be reused.